

Submittal Data Sheet

Project Information

Project _____ Approval _____


Specification

The NFPA 99 compliant digital, fully automatic manifold shall be a Tri-Tech Medical *Genesys™* series. No manual resetting of valves or levers shall be required. The unit shall switch from “Bank in Use” to “Reserve” bank without fluctuation in line delivery pressure. Simultaneously, the “Reserve in Use” alarm shall be triggered by the manifolds microprocessor. The manifold shall continue to provide gas, in the event of a power failure, until both banks are depleted. After the switchover, the “Reserve” bank shall then become the “Bank in Use”. The manifold microprocessor shall also trigger the “High Line Pressure” and “Low Line Pressure” alarms without the need for additional pressure switches or transducers. The manifold microprocessor shall also trigger the “Emergency Reserve in Use” and “Emergency Reserve Low” alarms when used with transducers supplied separately. The manifold shall be capable of being upgraded, to be used with high pressure (up to 3,000 psig) cylinders or for use at higher delivery pressures.

The control cabinet shall also incorporate economizer gas circuits for both banks. The economizer circuits will allow the head pressure of the reserve bank to be utilized instead of venting to atmosphere so long as there is sufficient system gas usage.

The control panel shall incorporate five large, green,” illuminated LED displays, for the Left Bank, the Right Bank, Delivery Pressure, Intermediate Pressure and Emergency Reserve Bank Pressure. Analog gauges are also provided so that all above pressure zones may be observed in the event of a power failure. The control panel shall also incorporate a set of LED’s for each bank, green for “Bank in Use”, amber for “Ready” and red for “Empty”.

All manifold regulators, piping and control switching equipment shall be cleaned for use with oxygen service and installed in a steel cabinet (weatherproof version available) to provide protection and minimize tampering






Features

- Input power 120 VAC, 50 to 60 Hz single point connection
- 72” flexible pigtailed between vaporizers and header
- 96” stainless double wall flexible pigtailed between cryogenic vessels and vaporizers.
- Economizer circuits for maximum efficiency of gas use.
- Unit of measure switching (psi, kPa, BAR).
- Dual line pressure regulators
- Expandable – additional portable bulk vessels may be added.
- Cabinet + headers weight 105 lbs.
- May be converted from low or medium pressure liquid portable bulk vessel use to use with high pressure cylinders.
- Line pressure sensor may be mounted inside the cabinet or remotely located to eliminate the need for a high/low pressure switch for master alarm operation.
- Flow capacity 3,500 scfh – (41.67 scfm) with a single portable bulk vessel.
- Packages available with 3,500 scfh non-continuous (12 hours on / 12 hours off) and continuous rated vaporizers.

System Requirements

Minimum inlet pressure	175 psig (pressure building circuits should be set @ 250 psig) vessels must have 350 psig relief valves
Physical space requirements	Continuous duty rated system – 16’ W x 8’ H x 45” D Non-Continuous duty rated system – 14’ W x 8’ H x 40” D



88-0967 Manifold Block CGA 440

Flow Information

Pressure at Inlet Block to Control Cabinet (just prior to changeover)	Delivery Pressure (allowing 5 psig drop)	Flow (SCFH / CFM)
95 (psig)	55 – 50	3,500 / 58.3

System Requirements

Minimum inlet pressure	175 psig (pressure building circuits should be set @ 250 psig) vessels must have 350 psig relief valves
Physical space requirements (includes vaporizers & 3 x 3 portable liquid vessels)	Continuous duty rated system – 16' W x 8' H x 45" D Non-Continuous duty rated system – 14' W x 8' H x 40" D
Cabinet Dimensions / Weight	26 ¼" H x 16 ¾" W x 9 ¼" D / 89 pounds Cabinet door is 28 ½" H x 19 ¼" W
Electrical Power Feed / Usage	A single point 120 volt AC 50 – 60 Hz / 45 W (0.4 amps)

How to Order

LU3522OX1H0000X	Weatherproof control cabinet with headers and two cryogenic transfer hose for 50 psi
LU3522OX2H0000X	Weatherproof control cabinet with headers and two cryogenic transfer hose for 80 psi
1X1CRYOKIT440	2 each 10' cryogenic liquid transfer hoses with fittings to connect to vaporizers (3/4 F NPT assumed)
2X2CRYOKIT440	4 each 10' cryogenic liquid transfer hoses with fittings to connect to vaporizers (3/4 F NPT assumed)
3X3CRYOKIT440	6 each 10' cryogenic liquid transfer hoses with fittings to connect to vaporizers (3/4 F NPT assumed)
88-0967	Manifold block – CGA 440
PSM-24	DISS Union/Gauge & Transducer Assembly - Oxygen
RWP-9-3S or V	Manifolds - High Pressure Reserve w/ Port – 3 cyl's (S= staggered V = Vertical Crossover)
RWP-9-4S or V	Manifolds - High Pressure Reserve w/ Port – 4 cyl's (S= staggered V = Vertical Crossover)
RWP-9-5S or V	Manifolds - High Pressure Reserve w/ Port – 5 cyl's (S= staggered V = Vertical Crossover)
RWP-9-6S or V	Manifolds - High Pressure Reserve w/ Port – 6 cyl's (S= staggered V = Vertical Crossover)
14-3001	2,500 psi transducer with 12' wire cable for Emergency Reserve Low Alarm
14-3002	500 psi transducer with 8' wire cable for Emergency Reserve in Use
CV-050F	Check valve 1/2" F NPT x 1/2 F NPT
600NCH	600 scfh per hour NON-CONTINUOUS duty rated vaporizer – hang mount TQ09-22979
1200 NC	1,200 scfh per hour NON-CONTINUOUS duty rated vaporizer – floor mount