

*Medical Gas Coupler Outlet  
To Fit  
Fairfield™ Rail Systems*

**Does your facility have leaking or broken rail outlets that look like these?**



**You can now replace or repair your outlets!**



FF222

FRF3605

FF124



Exclusive cartridge style primary check valve – for fast, easy repairs!

Gas Service	A Overall Body Length (inches)	B Back Body Thread Length (inches)	Complete Replacement Bodies Part #	Replacement Cartridge Part #	Cartridge Replacement Tool Part #
Oxygen	1.75	.375	FF124	F1024ARC	F9900CT
Oxygen	2.00	.275	FF224	F1024ARC	F9900CT
Oxygen	1.88	.493	Call for Assistance	F1024ARC	F9900CT
Vacuum	2.00	.375	FF122	F1022RC	F9900CT
Vacuum	2.00	.275	FF222	F1022RC	F9900CT
Vacuum	2.16	.493	FRF3605	F1022RC	F9900CT
Air	2.00	.375	FF116	F1000RC	F9900CT
Air	2.00	.275	FF216	F1000RC	F9900CT
Air	2.16	.493	FRF3606	F1000RC	F9900CT

*Contact Tri-Tech for assistance – we'll help you select the correct item!*